

Technologies for microscopic image analysis

显微图像分析技术

第5课

Takao, Daisuke 高尾大辅

dtakao@mail.hzau.edu.cn

Overview of this course

第一章 光学显微术 Optical microscopy

1. History of microscopy
2. Principles of optics
3. Bright field microscopy
4. Fluorescence microscopy

第二章 图像数据处理与分析 Image data processing and analysis

1. Processing of digital image data
2. Quantitative image analysis

3. Application of ImageJ/Fiji software

4. Image processing and analysis with the Python coding environment

5. Bioimage informatics



第三章 先进的显微术 Advanced microscopy

1. Advanced imaging technology: from single molecules to living animals
2. Super-resolution microscopy
3. Other bioimaging techniques

第四章 人工智能的应用 Application of AI

1. Overview of AI technology
2. AI for bioimaging

Image processing and analysis using Fiji/ImageJ



The screenshot shows the 'Results' window in Fiji/ImageJ, displaying a table of measurements for the 19 objects. The table has columns for 'Area', 'Mean', 'XM', and 'YM'. The measurements are as follows:

	Area	Mean	XM	YM
5	242.468	255	127.618	55.575
6	224.958	255	75.566	88.039
7	257.725	255	50.160	96.059
8	239.605	255	99.060	103.044
9	181.299	255	49.589	111.304
10	221.061	255	87.347	116.776
11	191.393	255	53.379	124.238
12	195.242	255	68.254	128.631
13	147.687	255	33.749	138.048
14	252.514	255	189.439	166.199
15	213.222	255	163.293	170.184
16	194.209	255	159.298	183.599
17	207.119	255	130.516	190.925
18	170.643	255	31.911	193.760
19	179.046	255	7.964	211.188

Tutorial: Image processing and analysis with Python


- Use image data from a published paper (Nagao et al., MBoC, 2020)
- Learn methods of segmentation, measurements, drawing graphs, etc

Preparation of the Python coding environment

Install and setup Anaconda (Python + other powerful packages for data science, informatics, and AI)

Not sure which one to use? – Anaconda is recommended for most of biologists

Python



Anaconda



► Home

▼ Anaconda Distribution

Installation

Installing on Windows

Installing on macOS

Installing on Linux

Installing on AWS Graviton2 (arm64)

Installing on Linux-s390x (IBM Z)

Installing on Linux POWER

Installing in silent mode

Installing for multiple users

Verifying your installation

Anaconda installer file hashes

Updating from older versions

Using Anaconda on older operating systems

Uninstalling Anaconda Distribution

User guide

Reference

End User License Agreement - Anaconda Distribution

► Anaconda Navigator



Installation

Review the system requirements listed below before installing Anaconda Distribution. If you don't want the hundreds of packages included with Anaconda, [install Miniconda](#), a mini version of Anaconda that includes just conda, its dependencies, and Python.



Tip

Looking for Python 3.5 or 3.6? See our [FAQ](#).

System requirements

- License: Free use and redistribution under the terms of the [EULA for Anaconda Distribution](#).
- Operating system: Windows 8 or newer, 64-bit macOS 10.13+, or Linux, including Ubuntu, RedHat, CentOS 7+, and others.
- If your operating system is older than what is currently supported, you can find older versions of the Anaconda installers in our [archive](#) that might work for you. See [Using Anaconda on older operating systems](#) for version recommendations.
- System architecture: Windows- 64-bit x86; MacOS- 64-bit x86 & M1; Linux- 64-bit x86, 64-bit aarch64 (AWS Graviton2), 64-bit Power8/Power9, s390x (Linux on IBM Z & LinuxONE).
- Minimum 5 GB disk space to download and install.

On Windows, macOS, and Linux, it is best to install Anaconda for the local user, which does not require administrator permissions and is the most robust type of installation. However, with administrator permissions, you can install Anaconda system wide.


- Installing on Windows
- Installing on macOS

- Installing on Linux
- Installing on AWS Graviton2 (arm64)
- Installing on Linux-s390x (IBM Z)
- Installing on Linux POWER
- Installing in silent mode
- Installing for multiple users

<https://docs.anaconda.com/anaconda/install/index.html>

Preparation of the Python coding environment

Concept of virtual environments



Package C only works with Python 3.6
and requires Package B v1.0

Multiple virtual environments can be constructed in one computer



Virtual environment 1: Cell segmentation-related packages installed

Virtual environment 2: Multivariate analysis-related packages installed

⋮

Preparation of the Python coding environment

1. Open Anaconda Navigator



Preparation of the Python coding environment

The screenshot shows the Anaconda Navigator desktop application. The interface includes a top menu bar with 'File' and 'Help', a sidebar on the left with navigation options, and a main workspace displaying various application tiles. Three blue callout boxes with white text and numbers are overlaid on the interface:

- 1** Points to the 'Home' button in the sidebar.
- 2** Points to the 'tutorial' dropdown menu in the top right of the main workspace.
- 3** Points to the 'JupyterLab' tile in the main workspace, which has a blue border and a 'Launch' button.


The sidebar on the left contains the following elements:

- Home** (with a house icon)
- Environments** (with a cube icon)
- Learning** (with a book icon)
- Community** (with a group of people icon)
- Anaconda Notebooks** (with a 'New!' badge and a 'Learn More' button)
- Documentation** (with a link icon)
- Anaconda Blog**
- Social media icons for Twitter, YouTube, and GitHub.

The main workspace displays a grid of application tiles under the 'tutorial' channel:

- DataSpell**: An IDE for exploratory data analysis and prototyping machine learning models. It combines the interactivity of Jupyter notebooks with the intelligent Python and R coding assistance of PyCharm in one user-friendly environment. **Install** button.
- JupyterLab 3.5.3**: An extensible environment for interactive and reproducible computing, based on the Jupyter Notebook and Architecture. **Launch** button.
- jupyter**: Web-based, interactive computing notebook environment. Edit and run human-readable docs while describing the data analysis. **Launch** button.
- VS Code 1.76.2**: Streamlined code editor with support for development operations like debugging, task running and version control.
- Datalore**: Kick-start your data science projects in seconds in a pre-configured environment. Enjoy coding assistance for Python, SQL, and... **Launch** button.
- IBM Watson Studio Cloud**: IBM Watson Studio Cloud provides you the tools to analyze and visualize data, to cleanse and shape data, to create and train machine... **Launch** button.

Preparation of the Python coding environment



Preparation of the Python coding environment

If you opened a file, shut down the kernel before closing the file

1

2. "Shut Down All" or close each kernel

```
[1]: import pandas as pd
import seaborn as sns
```

```
[2]: # set directory and file names
dir = "C:/Users/takao/Desktop" # change to your directory where the data file is
file1 = "data.csv" # import file name
file2 = "plot.pdf" # export file name

# import CSV file as dataframe
df = pd.read_csv(os.path.join(dir, file1)) # os.path.join is a convenient function to handle directories

# show imported dataframe
df
```

	0	125	250	500
0	66.499	115.198	158.120	125.026
1	55.148	265.979	298.091	1065.642
2	73.994	20.836	594.527	565.073
3	67.255	283.237	25.313	313.071
4	252.093	60.310	577.782	592.740
5	472.591	312.830	612.532	651.453
6	219.044	106.554	136.916	659.412
7	290.884	184.523	268.792	258.895
8	144.899	0.477	129.258	997.445
9	1514.207	5.116	510.290	671.624
10	601.792	73.494	1220.035	56.836
11	2.425	289.284	204.973	369.867

- For another option, you can install and setup VS Code (Visual Studio Code)
- You will still need to install Anaconda or Python first
- Setting up VS Code is a little complicated but VS Code will provide a more professional environment with a similar look


17 57.691 345.894 332.547 462.723

Simple 0 1 Python 3 (ipykernel) | Idle Mode: Command Ln 1, Col 1 Beeswarm-Boxplot_tutorial.ipynb

Data

- Hoechst_mean.csv
- HeLa_GM130_CENP-F (20 TIFF files of #011–030)

Jupyter notebook (ipynb) and related files

- HelloWorld.ipynb
 - Beeswarm-Boxplot_tutorial.ipynb
 - Image_Quantification_tutorial_1.ipynb
 - Image_Quantification_tutorial_2.ipynb
 - Image_Quantification_tutorial_3.ipynb
- 
- Proceed through the files in this order
- tutorial.yaml (used for setting up the environment; slide 8, step 5)

“hello world”

Write in a “cell” and run it by pressing “Shift” + “Enter/Return”

```
[1]: 1+1
[1]: 2

[2]: 1 + 1 # space will be ignored
[2]: 2

[3]: a = 1 # assign "1" to a variable "a"

[4]: print(a) # show the content of "a"
1

[5]: number = 2 # the names of variables are flexible
print(number) # show the content of "number"
2

[6]: b = a + number # variables can be used for calculations
print(b) # show the result
3

[7]: print("hello world") # show texts
hello world

[8]: text = "hello world" # texts can be assigned to variables
print(text)
hello world

[9]: number + text # this will result in error

-----
TypeError                                Traceback (most recent call last)
Cell In [9], line 1
----> 1 number + text

TypeError: unsupported operand type(s) for +: 'int' and 'str'

[10]: import numpy as np # installed packages can be imported like this


[11]: list1 = [0, 1, 2, 3]

[12]: np.mean(list1) # calculate the mean of list1 using the "mean" function of numpy
[12]: 1.5

[13]: # calculate the mean, min, and max values of list1 and show the results
list1_mean = np.mean(list1)
list1_min = np.min(list1)
list1_max = np.max(list1)
print("mean = ", list1_mean, "; min = ", list1_min, "; max = ", list1_max)

mean = 1.5 ; min = 0 ; max = 3
```

Plot measured data



	A	B	C	D
1	G1/S	G2		
2	121.248	135.7469		
3	112.9679	144.4894		
4	111.0115	152.8638		
5	105.0882	108.355		
6	93.8516	135.358		
7	95.29149	130.8957		
8	94.09235	149.2067		
9	94.04405	153.8039		
10	100.0526	186.7222		

Nagao et al., MBoC, 2020

HeLa cells

- Mean fluorescence intensities of Hoechst (DNA) were measured and listed in Hoechst_mean.csv.
- Data were categorized according to cell cycle phases (G1/S and G2)